tutorial

April 6, 2023

```
[27]: import sys
    sys.path.append('./test')
    import pylab as pl
    import numpy as np
    import xarray as xr
    from rhkpy import rhkpy
```

1 Loading dI/dV map

```
[28]: filename_map = './test/dI-dV_ABC_Graphite-Sample5_9K_2021_08_24_11_54_33_632.

sm4'

specmap = rhkpy.stmdata(filename_map)
```

There are two xarray Datasets in the file: image and spectra. One can access these using the dot . notation

```
[9]: specmap.image
[9]: <xarray.Dataset>
                      (x: 512, y: 512, scandir: 2)
     Dimensions:
     Coordinates:
       * x
                      (x) float64 -563.6 -563.4 -563.2 -563.0 ... -464.2 -464.0 -463.8
                      (y) float64 -828.9 -828.7 -828.5 -828.3 ... -729.5 -729.3 -729.1
       * y
       * scandir
                      (scandir) <U8 'forward' 'backward'
     Data variables:
                     (x, y, scandir) float64 0.2228 0.5587 0.3199 ... 1.126 1.268
         topography
         current
                      (x, y, scandir) float64 105.2 95.93 108.6 ... 99.23 92.91 96.35
         lia
                      (x, y, scandir) float64 -0.05732 -0.1006 ... -0.09234 -0.1453
     Attributes:
         filename:
                             ./test/dI-dV_ABC_Graphite-Sample5_9K_2021_08_24_11_54_...
         xoffset:
                             -563.5723326555265
                             -828.9023483917802
         yoffset:
         xoffset_units:
         yoffset_units:
                             nm
         bias:
                             0.5
         bias units:
         setpoint:
                             100.0000013351432
```

setpoint units: pA

measurement date: 08/23/21 measurement time: 19:46:50

Both the spectra and image have attributes, which can be accessed like so:

```
[39]: specmap.image.attrs
[39]: {'filename': './test/dI-dV_ABC_Graphite-Sample5_9K_2021_08_24_11_54_33_632.sm4',
       'xoffset': -563.5723326555265,
       'yoffset': -828.9023483917802,
       'xoffset_units': 'nm',
       'yoffset_units': 'nm',
       'bias': 0.5,
       'bias units': 'V',
       'setpoint': 100.0000013351432,
       'setpoint units': 'pA',
       'measurement date': '08/23/21',
       'measurement time': '19:46:50'}
     image has data variables: 'topography', 'current' and 'lia'. It has coordinates: 'x', 'y', 'scandir'
     Plotting the topography data, we select the forward scan direction. This can be done using the
     isel function
[40]: fwscan = specmap.image.isel(scandir=0)
     Now we can see that the coordinate scandir is gone, leaving only x and y
[41]: fwscan
[41]: <xarray.Dataset>
      Dimensions:
                       (x: 512, y: 512)
      Coordinates:
                       (x) float64 -563.6 -563.4 -563.2 -563.0 ... -464.2 -464.0 -463.8
        * x
                       (y) float64 -828.9 -828.7 -828.5 -828.3 ... -729.5 -729.3 -729.1
        * y
          scandir
                       <U8 'forward'
      Data variables:
                       (x, y) float64 0.2228 0.3199 0.4098 0.348 ... 1.181 1.199 1.126
          topography
                       (x, y) float64 105.2 108.6 94.26 95.62 ... 106.6 98.58 92.91
          current
                       (x, y) float64 -0.05732 -0.09757 0.2002 ... 0.01556 -0.09234
          lia
      Attributes:
          filename:
                               ./test/dI-dV_ABC_Graphite-Sample5_9K_2021_08_24_11_54_...
          xoffset:
                               -563.5723326555265
          yoffset:
                               -828.9023483917802
          xoffset_units:
                               nm
          yoffset_units:
                               nm
          bias:
                               0.5
          bias units:
                               V
```

```
setpoint: 100.0000013351432
setpoint units: pA
measurement date: 08/23/21
measurement time: 19:46:50
```

Selecting the 'topography' DataArray from the set. This has coordinates 'x' and 'y' and values as a numpy array, which can be accessed by fwscan['topography'].data

```
[18]: fwscan['topography']
[18]: <xarray.DataArray 'topography' (x: 512, y: 512)>
      array([[ 0.22277586, 0.31985992, 0.40975636, ..., 0.1738827,
              0.25653668, 0.21637588],
             [0.42118773, 0.49016015, 0.49237385, ..., 0.55386099,
              0.4927506 , 0.47112197],
             [0.66517782, 0.67569847, 0.68743548, ..., 0.68594219,
              0.71097722, 0.65945614
             [-0.04212871, -0.02539763, 0.02653675, ..., 1.07600668,
              1.05096369, 1.07442425],
             [-0.1272415, -0.10926353, -0.07156714, ..., 1.10464527,
              1.05920205, 1.02997055],
             [0.05461745, 0.07771955, 0.08524178, ..., 1.18111791,
              1.19945448, 1.12606963]])
      Coordinates:
                   (x) float64 -563.6 -563.4 -563.2 -563.0 ... -464.2 -464.0 -463.8
        * x
                   (y) float64 -828.9 -828.7 -828.5 -828.3 ... -729.5 -729.3 -729.1
         scandir <U8 'forward'
      Attributes:
         units:
                      nm
         long units: nanometer
[43]: fwscan['topography'].data
[43]: array([[ 0.22277586, 0.31985992, 0.40975636, ..., 0.1738827,
              0.25653668, 0.21637588],
             [0.42118773, 0.49016015, 0.49237385, ..., 0.55386099,
              0.4927506 , 0.47112197],
             [0.66517782, 0.67569847, 0.68743548, ..., 0.68594219,
              0.71097722, 0.65945614],
             [-0.04212871, -0.02539763, 0.02653675, ..., 1.07600668,
               1.05096369, 1.07442425],
             [-0.1272415, -0.10926353, -0.07156714, ..., 1.10464527,
              1.05920205, 1.02997055],
             [0.05461745, 0.07771955, 0.08524178, ..., 1.18111791,
               1.19945448, 1.12606963]])
```

The coordinates in the Dataset are accessed by .coords, for example in the 'x' direction:

fwscan.coords['x']. The values can be directly acessed as: fwscan.coords['x'].data

[46]: fwscan.coords['x']

[46]: <xarray.DataArray 'x' (x: 512)>
 array([-563.572333, -563.37702 , -563.181708, ..., -464.158269, -463.962956,
 -463.767644])

Coordinates:

* x (x) float64 -563.6 -563.4 -563.2 -563.0 ... -464.2 -464.0 -463.8 scandir <U8 'forward'

Attributes:

units: nm


long units: nanometer

Plotting the 'topography' DataArray. When plotting the 'topography' the data has been plane and line fitted.

The coordinates x and y are in the absolute tip positions

[16]: fwscan['topography'].plot()

[16]: <matplotlib.collections.QuadMesh at 0x22cdf402c20>



Let's see the spectra Dataset


```
[29]: specmap.spectra
[29]: <xarray.Dataset>
     Dimensions:
                        (bias: 501, specpos_x: 32, specpos_y: 32, repetitions: 1,
                        biasscandir: 2)
      Coordinates:
        * bias
                        (bias) float64 0.5 0.498 0.496 0.494 ... -0.496 -0.498 -0.5
        * specpos_x
                        (specpos_x) float64 -465.3 -468.5 -471.6 ... -559.1 -562.2
                        (specpos_y) float64 -730.7 -733.8 -736.9 ... -824.4 -827.5
        * specpos_y
                       (repetitions) int32 0
        * repetitions
        * biasscandir
                       (biasscandir) <U5 'left' 'right'
     Data variables:
                        (bias, specpos_x, specpos_y, repetitions, biasscandir) float64
          lia
                        (bias, specpos_x, specpos_y, repetitions, biasscandir) float64
          current
                        (specpos_x, specpos_y) float64 -465.3 -468.5 ... -559.1 -562.2
          х
                        (specpos_x, specpos_y) float64 -730.7 -730.7 ... -827.5 -827.5
      Attributes:
          filename:
                              ./test/dI-dV ABC Graphite-Sample5 9K 2021 08 24 11 54 ...
          bias:
                              0.5
                              V
          bias units:
                              100.0000013351432
          setpoint:
          setpoint units:
          measurement date:
                             08/23/21
          measurement time:
                             19:50:01
```

Plotting the dI/dV spectra.

We can use the sel function for this, by specifying the coordinates in nm. Finally selecting the 'lia' data variable to plot as a function of the bias coordinate. To do this we first select the one value in repetitions coordinate, because there are only two spectra (left and right bias scan) at each tip position. We also can use the isel and sel functions together.


```
[36]: specmap2 = specmap.spectra.isel(repetitions=0) specmap2.isel(biasscandir=0).sel(specpos_x=-782, specpos_y=-521,u omethod='nearest')['lia'].plot()
```

[36]: [<matplotlib.lines.Line2D at 0x22ce4e5f7c0>]



Plot both the current and Lock-In






Plot the dI/dV value at zero bias.

```
[33]: specmap.spectra.isel(repetitions=0, biasscandir=0).sel(bias=0, ⊔

→method='nearest')['lia'].plot()
```


[33]: <matplotlib.collections.QuadMesh at 0x22cdf2046d0>



Average over the forward and backward directions of the bias sweep

```
[35]: specmap.spectra.mean(dim=['repetitions', 'biasscandir']).sel(bias=0, use thod='nearest')['lia'].plot()
```

[35]: <matplotlib.collections.QuadMesh at 0x22ce6b5bcd0>



Let's do something more fancy, like plotting spectra across a line, parallel to the 'y' axis at 'x'=-542 nm

```
[48]: specmap_avg = specmap.spectra.mean(dim=['repetitions', 'biasscandir'])
```

The selected data now has coordinates of bias and specpos_y.

[51]: spec_along_line

[51]: <xarray.Dataset>

Dimensions: (bias: 501, specpos_y: 32)

Coordinates:

- * bias (bias) float64 0.5 0.498 0.496 0.494 ... -0.496 -0.498 -0.5 specpos_x float64 -543.5
- * specpos_y (specpos_y) float64 -730.7 -733.8 -736.9 ... -821.3 -824.4 -827.5 Data variables:


lia (bias, specpos_y) float64 18.12 19.5 15.82 ... 17.33 16.82 15.61 current (bias, specpos_y) float64 122.6 130.2 95.7 ... -136.6 -129.6 x (specpos_y) float64 -465.3 -468.5 -471.6 ... -556.0 -559.1 -562.2

y (specpos_y) float64 -808.8 -808.8 -808.8 ... -808.8 -808.8 -808.8

We can plot this on a density plot.

[53]: spec_along_line['lia'].plot()

[53]: <matplotlib.collections.QuadMesh at 0x22ce7177c70>



[]: